

# Lab::VISA

## & friends

slides © Daniel Schröder and Andreas K. Hüttel

# VISA, GPIB, etc.

- Instruments can be connected in various ways: Serial port, GPIB, VXI, TCP/IP, etc.
- **GPIB** (hardware and software)
  - GPIB (IEEE488): Standard by Hewlett-Packard
  - Physical layer IEEE488.1
  - Command layer IEEE488.2
  - SCPI (Standard Commands for Programmable Instruments)
- **VISA** (software)
  - Virtual Instrument Software Architecture
  - VXI, GPIB, serial, or computer-based instruments
  - NI-VISA library is one implementation of the VISA standard

# Lab::VISA

- Lab::VISA is software to control instruments in the lab via VISA (e.g. voltage sources, multimeters, etc.)
- Alternative to LabView, Gpplus, etc.
- But very different: no GUI, just API

# Lab::VISA design goals

- **Flexible**
  - Allow any kind of measurement procedure
  - Control anything that has GPIB or serial connection
- **Safe**
  - Make sure you never drive voltage sources too fast and destroy your gates
- **Helpful**
  - Automatically collect additional information about the data (“metadata”)
  - User should have to type in additional information only once and then never again

# Lab::VISA is in Perl

- Lab::VISA is written in Perl and C
- To be used from programs written in Perl
- Perl: interpreted scripting language
- Runs on almost every OS
- Extremely good in reading data files, manipulating data, etc.
- Allows to write quick and dirty scripts that get the job done
- Also allows to write clean fullsize programs


=> ideal for experimental physics

# Lab::VISA architecture

- Divided into three parts
- Build on top of each other
- Provide increasing comfort
- Measurement scripts can be based on any of these stages
- **Lab::VISA** is the lowest layer. It makes the NI-VISA library accessible from perl and therefor allows to make any standard VISA call
- **Lab::Instrument** package makes communication with instruments easier by silently handling the protocol involved
- **Lab::Tools** is the highest abstraction layer. Provides support for writing good measurement scripts. Offers means of saving data and related meta information to disk, plotting data, etc.

# Lab::VISA architecture

```
your (tiny) measurement program / script
```



# Example 1: Use only Lab::VISA

- looks like example programs from the NI VISA manual
  - protocol overhead is pain in the neck. **NOT GOOD.**

```
#!/usr/bin/perl

use strict;
use Lab::VISA;

# Initialize VISA system and
# Open default resource manager
my ($status,$default_rm)=Lab::VISA::viOpenDefaultRM();
if ($status != $Lab::VISA::VI_SUCCESS) {
    die "Cannot open resource manager: $status";
}

# Open one resource (an instrument)
my $gpib=21;      # we want to open the instrument
my $board=0;      # with GPIB address 21
                  # connected to GPIB board 0 in our computer

my $resource_name=sprintf("GPIB%u::%u::INSTR",$board,$gpib);

($status, my $instr)=Lab::VISA::viOpen(
    $default_rm,    # the resource manager session
    $resource_name, # a string describing the
    $Lab::VISA::VI_NULL,# access mode (no special mode)
    $Lab::VISA::VI_NULL # time out for open (no time out)
);

if ($status != $Lab::VISA::VI_SUCCESS) {
    die "Cannot open instrument $resource_name. status: $status";
}

(...)
```

```
(...)

# Now we are going to send one command and read the result.

# We send the simple SCPI command "*IDN?" which asks the instrument
# to identify itself. Of course the instrument must support this
# command, in order to make this example work.
my $cmd="*IDN?";
($status, my $write_cnt)=Lab::VISA::viWrite(
    $instr,      # the session identifier
    $cmd,        # the command to send
    length($cmd) # the length of the command in bytes
);
if ($status != $Lab::VISA::VI_SUCCESS) {
    die "Error while writing: $status";
}

# Now we will read the instruments reply
($status,      # indicates if the operation was successful
my $result,    # the answer string
my $read_cnt)= # the length of the answer in bytes
    Lab::VISA::viRead(
        $instr, # the session identifier
        300     # read 300 bytes
    );
if ($status != $Lab::VISA::VI_SUCCESS) {
    die "Error while reading: $status";
}
# The result string will be 300 bytes long, but only $read_cnt
# bytes are part of the answer. We cut away the rest.
$result=substr($result,0,$read_cnt);

print $result;
```



# Example 2: Use Lab::Instrument

- Much nicer!

```
#!/usr/bin/perl
```

```
use strict;
```

```
use Lab::Instrument::HP34401A;
```

```
my $gpib=21;          # we want to open the instrument
```

```
my $board=0;          # with GPIB address 21
```

```
                    # connected to GPIB board 0 in our computer
```

```
# Create an instrument object
```

```
my $hp=new Lab::Instrument::HP34401A($board,$gpib);
```

```
# Use the id method to query the instruments ID string
```

```
my $result=$hp->id();
```

```
print $result;
```

# Example 3: My first real measurement!

- How about doing a gate sweep?
- That's easy!

```
#!/usr/bin/perl

use strict;
use Lab::Instrument::HP34401A;
use Lab::Instrument::Yokogawa7651;

# Connect to voltage source
my $yoko=new Lab::Instrument::Yokogawa7651({
    GPIB_address => 14,
    gate_protect => 0,
});

# Connect to multimeter
my $hp=new Lab::Instrument::HP34401A(0,21);

# Sweep Yokogawa from 0 to 1 volt (10 steps)

for (my $volt=0; $volt<=1; $volt=$volt+0.1) {
    # set Yokogawa
    $yoko->set_voltage($volt);

    # wait a second
    sleep(1);


    # read multimeter
    my $vmeas=$hp->read_voltage_dc();

    # print values
    print $volt,"\t",$vmeas,"\n";
}
```

# Lab::Tools

- Provides additional tools to write better measurement scripts
- Store **metadata** alongside **data**
  - date and time
  - settings of additional instruments
  - ratio of voltage divider
  - color of the shirt you are wearing
  - everything that might be important for a later interpretation of the data
- Don't repeat yourself
  - Use above collected information automatically
  - Automatically plot data with correct axes, scaling, labels etc.

# Metadata philosophy



## Constants

- $R = 1000 \text{ Ohm}$

## Columns

- $C1: V_{\text{bias}}$
- $C2: V_{\text{resistor}}$
- unit
- description...

1.5	+8.19300500E-01
1.4	+7.23413000E-01
1.3	+6.28083900E-01
1.5	+8.19300500E-01

## Axes

- Axis “bias voltage”  
( $C1$ )
- Axis “diode current”  
( $C2 / R$ )
- Axis “diode resistance”  
( $R * (C1 / C2 - 1)$ )
- unit
- expression
- description...

## Plots

- Plot “diode current”  
axis “diode current” over  
axis “bias voltage”
- Plot “diode resistance”  
axis “diode resistance” over  
axis “bias voltage”
- logscale
- grid
- ranges...

# Example 4: Using Lab::Measurement

- metadata, live plot

```
#!/usr/bin/perl

use strict;
use Lab::Instrument::KnickS252;
use Lab::Instrument::HP34401A;
use Lab::Measurement;
use Time::HiRes qw/usleep/;

#####

my $start_voltage = 1.5;
my $end_voltage   = -3.5;
my $step          = -0.05;

my $knick_gpib    = 14;
my $hp_gpib       = 21;

my $sample        = "Zenerdiode";
my $title         = "Messung mit Lab::Measurement";
my $comment       = <<COMMENT;
Reihenschaltung aus Widerstand 1kOhm und Zenerdiode.
COMMENT

#####

my $knick=new Lab::Instrument::KnickS252({
    'GPIB_board' => 0,
    'GPIB_address' => $knick_gpib,
    'gate_protect' => 0,
});

my $hp=new Lab::Instrument::HP34401A(0,$hp_gpib);

my $measurement=new Lab::Measurement(
    sample => $sample,
    title => $title,
    filename_base => 'zener_kennlinie',
    description => $comment,

    live_plot => 'diode current',
```

```
constants => [
    {
        'name' => 'R',
        'value' => '1000',
    },
],
columns => [
    {
        'unit' => 'V',
        'label' => 'V_{bias}',
        'description' => 'Bias Voltage',
    },
    {
        'unit' => 'V',
        'label' => 'Amplifier output',
        'description' => 'Voltage drop on serial resistor',
    },
],
axes => [
    {
        'unit' => 'V',
        'expression' => '$C0',
        'label' => 'V_{bias}',
        'description' => 'Bias voltage',
        'min' => ($start_voltage < $end_voltage)
            ? $start_voltage
            : $end_voltage,
        'max' => ($start_voltage < $end_voltage)
            ? $end_voltage
            : $start_voltage,
    },
    {
        'unit' => 'mA',
        'expression' => '1000*($C1/R)',
        'label' => 'I_{diode}',
        'description' => 'Current through diode',
        'min' => '-1.2',
        'max' => '1.2',
    },
],
```

```
{
    'unit' => 'Ohm',
    'expression' => 'R * ($C0/$C1 - 1)',
    'label' => 'R_{diode}',
    'description' => 'Diode resistance',
},
],
plots => {
    'diode current' => {
        'type' => 'line',
        'xaxis' => 0,
        'yaxis' => 1,
        'grid' => 'xtics ytics',
    },
    'diode resistance' => {
        'type' => 'line',
        'xaxis' => 0,
        'yaxis' => 2,
        'grid' => 'xtics ytics',
        'logscale' => 'y',
    },
},
);

$measurement->start_block();

for (
    my $volt = $start_voltage;
    ($volt - $end_voltage) / $step < 0.5;
    $volt += $step
){
    $knick->set_voltage($volt);
    usleep(500000);
    my $meas = $hp->read_voltage_dc(10,0.0001);
    $measurement->log_line($volt,$meas);
}

my $meta = $measurement->finish_measurement();
```

# More utilities: `plotter.pl`, `make_overview.pl`

- **`plotter.pl`:**
  - Reads measurement file, list available plots (axes etc)
  - Creates postscript or pdf output
- **`make_overview.pl`:**
  - Reads all measurements in a directory
  - Generates a postscript or pdf file with plots for each measurement, including the metadata (i.e. constants, parameters, color of your shirt) (via LaTeX)
  - Great for completing your lab book

# Other Lab::VISA features

- “**Gate protect**” safety mechanism
  - Makes sure that no voltage is changed too fast
  - Big voltage steps are automatically split into small, slow steps
- **Date/time** handling
  - Date/time column can contain timestamp for every data point
  - Plot data as function of time
- Measurements with **higher dimensionality**
  - Each trace/sweep/line is a “block”
  - Two-dimensional plots, selections of traces, etc.